

https://www.innovationforever.com

Journal of

Modern Nursing Practice and Research

ISSN 2708-2202 (Online)

Open Access

Short Communication

Feasibility of Non-contact Noise Cancellation for Medical Device Alarms in Intensive Care Unit: An *In Vitro* Experiment

Nao Sonoda^{1*}, Kiyotaka Fujii², Akiko Morimoto¹

¹Graduate School of Nursing, Osaka Metropolitan University, Osaka, Japan

²School of Allied Health Sciences, Kitasato University, Kanagawa, Japan

***Correspondence to: Nao Sonoda, PhD, RN,** Graduate School of Nursing, Osaka Metropolitan University, 3-7-30 Habikino, Habikino-city, Osaka, 583-8555, Japan; Email: sonoda@omu.ac.jp

Received: June 4, 2024 Revised: July 11, 2024 Accepted: July 23, 2024 Published: August 6, 2024

Abstract

Objective: Medical device alarms in intensive care units (ICUs) are essential but their constant noise contributes to sleep disturbance and delirium in patients, jeopardizing their health and recovery. Therefore, improving the sound environment for patients in the ICU is crucial. This fundamental experiment investigated the potential of non-contact noise cancellation to reduce the noise level from medical device alarms in patients' ears. We tested whether non-contact noise cancellation using a loudspeaker reduced the noise level of medical device alarms.

Methods: In Experiment 1, we measured and analyzed the electrocardiogram sounds (volume 2, default setting) from a bedside monitor and the alarm sounds (volume 2, default setting) from a syringe pump. Negative-phase sound waves were subsequently generated for each alarm type. In Experiment 2, we tested whether non-contact noise cancellation using a loudspeaker reduced the noise level of each alarm.

Results: The electrocardiogram sound from the bedside monitor was highest around 1kHz, while, the alarm sound from the syringe pump peaked around 4kHz. Non-contact noise cancellation using the loudspeaker reduced the electrocardiogram sounds emitted from the bedside monitor by more than 15 decibel (dB) and the alarm sounds from the syringe pump by more than 10dB.

Conclusion: A 10dB change doubles or halves the perceived loudness, making the over 10dB reduction achieved for medical device alarms in this study highly significant. The results of this experimental study demonstrate the promising potential of noise-cancelling loudspeakers for improving patient comfort and well-being in ICUs.

Keywords: Noise, sound environment, non-contact noise cancellation, medical device alarm, ICU

Citation: Sonoda N, Fujii K, Morimoto A. Feasibility of Non-contact Noise Cancellation for Medical Device Alarms in Intensive Care Unit: An *In Vitro* Experiment. *J Mod Nurs Pract Res*, 2024; 4(3): 15. DOI: 10.53964/jmnpr.2024015.

1 INTRODUCTION

Medical device alarms in intensive care units (ICUs) are essential. However, their constant noise contributes to sleep disturbance and delirium in patients, jeopardizing their health and recovery^[1-6]. Patients with delirium are at an increased risk of dying during admission, longer stays in the hospital, and cognitive impairment after discharge^[7]. Therefore, it is important to improve the sound environment for patients in the ICU.

The World Health Organization recommends that noise levels (LAeq) should not exceed 35 decibel (dB) during the day and 30-40dB at night in hospital environments^[8]. However, LAeq in the ICU are frequently reported to be in the range of 45-70dB, with peaks exceeding 85dB^[4-6,9-13]. Although the alarm volume can be reduced from the default setting in the ICU, where patients are usually in critical condition, low-volume alarms may miss vital alerts, making medical professionals such as nurses anxious^[14]. Therefore, it is necessary to reduce LAeq perceived by patients without compromising medical professionals' awareness.

Noise reduction strategies include passive noise cancellation (e.g., earplugs), which physically prevents outside sounds from entering the ear, and active noise cancellation, which reduces the noise level in a particular area by sending a negative-phase sound. Earplugs are the commonest measure; however, the majority of studies related to the use of earplugs have found no benefit on objective sleep outcomes^[15]. Conversely, regarding active noise cancellation, an experimental study reported that commercially available noise-cancelling headphones reduced LAeq in an ICU environment by 6.8dB^[16]. However, noise-cancelling headphones are not suitable for patients in the ICU with poor general condition, who are at risk of developing pressure ulcers, or for those with head trauma.

Sound is a pressure wave consisting of alternating periods of compression and rarefaction. The active noise cancellation concept relies on the superposition principle: a noise-cancelling loudspeaker emits a sound wave with the same amplitude but with an inverted phase (negative-phase sound) relative to the original sound (positive-phase sound). The waves combine to form a new wave in a process called interference and effectively cancel each other out (Figure 1). A loudspeaker may enable non-contact noise cancellation in patients in the ICU.

As a fundamental experiment for developing non-contact noise cancellation to reduce the medical device alarm LAeq around patients' ears, our study tested whether non-contact noise cancellation using a loudspeaker reduced the noise level of medical device alarms.

2 METHODS

We measured and analyzed medical device alarm

sounds in Experiment 1 and created negative-phase sounds. In Experiment 2, we tested whether non-contact noise cancellation using a loudspeaker reduced the noise level of medical device alarms.

2.1 Experiment 1 (Characterization of Alarm Sounds) 2.1.1 Experiment Location

Experiment 1 was conducted in a simulated ICU environment (ICU Room at the Kobe International Living Medical Center, Hyogo, Japan) (Figure 2) in March 2020.

2.1.2 Medical Device Alarms

We measured the electrocardiogram sounds (volume 2, default setting) from a bedside monitor and alarm sounds (volume 2, default setting) from a syringe pump, representing high-frequency and simple alarm sounds commonly used in ICUs.

2.1.3 Measurement and Analysis

On the day of the measurement, no one was allowed on the floor except for the investigators. Additionally, all air conditioners and devices on the floor were turned off except for the medical devices to be measured. Measurements were recorded 100cm from the loudspeaker of each medical device (pillow position on the bed) (Figure 3). Further, each alarm sound was recorded for 10s using a sound level meter (NL-22, RION Corporation, Tokyo, Japan), and the LAeq (dB) and LAeq (dB) per frequency band were calculated. LAeq is the A-weighted equivalent continuous sound level (dB), and is a measure averaging noise exposure over a given period of time. Measurements were repeated thrice in the same environment.

2.2 Experiment 2 (Non-contact Noise Cancellation Test) 2.2.1 Experiment Location

Experiment 2 was conducted in an anechoic chamber (Hyogo, Japan) in December 2020 to experimentally and accurately evaluate noise reduction.


2.2.2 Positive-phase and Negative-phase Sounds

We used the electrocardiogram sounds from the bedside monitor and alarm sounds from the syringe pump obtained in Experiment 1. Negative-phase sounds were created for each alarm sound obtained in Experiment 1 using Audacity (version 2.1.3) (Audacity, New Delhi, India).

2.2.3 Measurement and Analysis

We tested whether non-contact noise cancellation using a loudspeaker reduced the noise level of each alarm. The following devices and software were used in the experiment: two loudspeakers (MSP3, Yamaha Corporation, Hamamatsu, Japan), a sound playback software (Live 9.7.7, Ableton Corporation, Berlin, Germany), a measurement microphone (Type4007, Bruel & Kjaer Corporation, Nærum, Denmark), and a sound measurement software (Smaart 8.4.3.1, Rational-Acoustics Corporation, Woodstock, CT, USA). The two

https://doi.org/10.53964/jmnpr.2024015



Loudspeaker

Figure 1. A noise reduction method using a noise-canceling loudspeaker.




Figure 3. Recording environment.

Table 1. Electrocardiogram Sounds From the BedsideMonitor and Alarm Sounds From the Syringe Pump

	Electrocardiogram Sounds (Volume 2, Default Setting) from the Dedside Monitor	Alarm Sounds (Volume 2, Default Setting) from the Syringe Pump	
LAeq for 10s (dB)			
1 st measurement	52.2	58.0	
2 nd measurement	52.1	58.9	
3 rd measurement	51.9	58.1	
Mean LAeq of three measurements $(dB)^{*}$	52.1±0.2	58.3±0.5	

Notes:*Mean ± standard deviation

loudspeakers were placed in parallel (Figure 4), with the measurement microphone in front of them. While one loudspeaker played positive-phase sounds (alarm sounds obtained in Experiment 1), the other loudspeaker played negative-phase sounds.

First, the positive-phase sound pressure level was set to 60dB at the 50cm point, and the LAeq was measured every 10cm between 50cm and 100cm points when only the positive-phase sound was played for 10s. Furthermore, the LAeq was measured at the same points with both positive-and negative-phase sounds playing simultaneously for 10s.



Figure 2. A simulated ICU environment.




Figure 4. Placement of two loudspeakers in the anechoic chamber.

Finally, the difference (dB) between "the LAeq when only the positive-phase sound was played" and "the LAeq when both positive- and negative-phase sounds were played simultaneously" was calculated.

3 RESULTS

3.1 Characterization of Alarm Sounds (Experiment 1)


Mean LAeq of three measurements for both electrocardiogram sounds and alarm sounds at 100cm from their respective loudspeakers exceeded 50dB (Table 1). Mean LAeq of three measurements for each frequency band is shown in Figure 5. The electrocardiogram sound (volume 2, default setting) from the bedside monitor was high in dB around 1kHz. Meanwhile, the alarm sound (volume 2, default setting) from the syringe pump was high in dB around 4kHz.

LAeq is the A-weighted equivalent continuous sound level (dB), and is a measure averaging noise exposure over a given period of time.

Measurements were performed at a distance of 100cm from the loudspeaker of each medical device (pillow position on the bed).

3.2 Non-contact Noise Cancellation Test (Experiment 2)

The degree of reduction in the noise level of each alarm due to non-contact noise cancellation using a loudspeaker



is shown in Table 2. The difference at the 100cm point was -16.2dB for the electrocardiogram sounds from the bedside monitor and -11.2dB for the alarm sounds from the syringe pump. Notably, both sounds experienced reductions exceeding 15dB and 10dB, respectively, at all points from 50cm to 100cm (average reductions: -17.3dB for electrocardiogram, -14.3dB for alarm).

LAeq is the A-weighted equivalent continuous sound level (dB), and is a measure averaging noise exposure over a given period of time.

4 DISCUSSION

The electrocardiogram sounds (volume 2, default setting) emitted from the bedside monitor were characterized by a high dB level approximately 1kHz, exceeding 50dB at a distance of 100cm from the monitor's loudspeaker. Notably, these sounds were reduced by more than 15dB owing to non-contact noise cancellation while using the loudspeaker. Similarly, the alarm sounds (volume 2, default setting) from the syringe pump, characterized by a high dB level approximately 4kHz and exceeding 50dB at a distance of 100cm from the loudspeaker of the syringe pump, were attenuated by over 10dB through non-contact noise cancellation.

Despite the increased emphasis on noise reduction in hospitals, hospital noise, particularly in ICUs, is a steadily worsening problem^[17]. Previous studies have reported that the main noise sources in the ICU are medical device alarms and the voices of medical professionals^[14,18-22]. Additionally, electronic sounds, such as alarms, have a greater negative effect on sleep disruption compared to human voices^[23]. Moreover, frequent sound peaks (abrupt sound elevation exceeding 10dB) are the commonest cause of noise-induced arousal^[24]. Notably, a 10dB change doubles or halves the perceived loudness, making the over 10dB reduction achieved for medical device alarms in this study highly significant. Our findings highlight the promising potential of noise-canceling loudspeakers to address this critical issue, and non-contact noise cancellation also has advantages in infection prevention.

	50cm	60cm	70cm	80cm	90cm	100cm
Electrocardiogram sounds from the bedside monitor						
LAeq for 10s when only the positive- phase sound was played (dB)	60.0	60.4	59.9	59.1	57.6	54.6
LAeq for 10s when positive- and negative-phase sounds were played simultaneously (dB)	40.1	40.1	43.7	43.7	41.7	38.4
Difference (dB)	-19.9	-20.3	-16.2	-15.4	-15.9	-16.2
Alarm sounds from the syringe pump						
LAeq for 10s when only the positive- phase sound was played (dB)	60.0	57.1	56.9	55.9	55.4	54.3
LAeq for 10s when positive- and negative-phase sounds were played simultaneously (dB)	42.5	41.5	37.3	45.6	43.7	43.1
Difference (dB)	-17.5	-15.6	-19.6	-10.3	-11.7	-11.2

Table 2. Degree of Reduction in the Noise Level of Each Alarm Due to Non-contact Noise CancellationUsing the Loudspeaker




Figure 6. Conceptualization of a future system that uses directional loudspeakers to reduce the noise level perceived by patients without inhibiting alerts to medical professionals.

The noise reduction recorded in this study surpassed the noise reduction (6.8dB reduction) achieved with commercially available noise-cancelling headphones in a previous experimental study^[16]. We attribute this improvement to the likely creation of negative-phase sounds specifically for each alarm sound. Moreover, in this study, the electrocardiogram sounds of the bedside monitor showed more noise reduction than the alarm sounds from the syringe pump. The inherent characteristics of active noise cancellation, which performs better at lower frequencies. This may explain why the alarm sounds from the syringe pump, characterized by high dB levels (approximately 4kHz), are less likely to be reduced.

Conventional loudspeakers diffuse sound in all directions. However, in recent years, directional loudspeakers, which transmit sound in a straight line, have been utilized in everyday situations. In the future, the use of directional loudspeakers sending sound to a precise location may create a system that reduces the noise level perceived by the patient without inhibiting alerts to medical professionals (Figure 6).

4.1 Limitations

This study had some limitations. It was a fundamental experiment for developing non-contact noise cancellation to reduce the noise level around patients' ears caused by medical device alarms. However, it is impossible to infer noise reduction in an actual ICU. Additionally, the alarm sounds tested in this study were limited.

5 CONCLUSION

In conclusion, this experimental study provides compelling evidence that non-contact noise cancellation using the loudspeaker reduced the electrocardiogram sounds

https://doi.org/10.53964/jmnpr.2024015

emitted from the bedside monitor by more than 15dB and the alarm sounds from the syringe pump by more than 10dB. The results of this study demonstrate the promising potential of noise-cancelling loudspeakers for improving patient comfort and well-being in ICUs. Further in vitro experiments are needed for other medical device alarms such as ventilators. Additionally, Further experiments in the actual ICU are needed to develop non-contact noise cancellation.

Acknowledgements

We were grateful to all investigators in this experimental study. We thanked Yuionn Corporation and Japan Noise Survey Company for their technical support. This work was supported by the Fund from the Terumo Life Science Foundation.

Conflicts of Interest

The authors declared no conflict of interest.

Ethical Statement

This study did not involve human participants and did not require approval from the institutional review board.

Author Contribution

Sonoda N contributed to conceptualization, data collection, analysis, funding acquisition, project administration, and writing of the original draft. Fujii K contributed to conceptualization, data collection, analysis, funding acquisition, review, and editing. Morimoto A contributed to conceptualization, data collection, funding acquisition, project administration, review, and editing. All authors thoroughly reviewed and endorsed the final version of the manuscript.

Abbreviation List

dB, Decibel ICU, Intensive care unit LAeq, Noise levels

References

- Xie H, Kang J, Mills GH. Clinical review: The impact of noise on patients' sleep and the effectiveness of noise reduction strategies in intensive care units. *Crit Care*, 2009; 13: 208.[DOI]
- [2] Figueroa-Ramos MI, Arroyo-Novoa CM, Lee KA et al. Sleep and delirium in ICU patients: a review of mechanisms and manifestations. *Intensive Care Med*, 2009; 35: 781-795.[DOI]
- [3] Sangari A, Emhardt EA, Salas B et al. Delirium Variability is Influenced by the Sound Environment (DEVISE Study): how changes in the intensive care unit soundscape affect delirium incidence. *J Med Syst*, 2021; 45: 76.[DOI]
- [4] Aydın Sayılan A, Kulakaç N, Sayılan S. The effects of noise levels on pain, anxiety, and sleep in patients. *Nurs Crit Care*, 2021; 26: 79-85.[DOI]

- [5] Bani Younis M, Hayajneh F, Alshraideh JA. Effect of noise and light levels on sleep of intensive care unit patients. *Nurs Crit Care*, 2021; 26: 73-78.[DOI]
- [6] Guisasola-Rabes M, Solà-Enriquez B, Vélez-Pereira AM et al. Noise levels and sleep in a surgical ICU. *J Clin Med*, 2022; 11: 2328.[DOI]
- [7] Salluh JI, Wang H, Schneider EB et al. Outcome of delirium in critically ill patients: systematic review and meta-analysis. *BMJ*, 2015; 350: h2538.[DOI]
- [8] Berglund B, Lindvall T, Schwela DH. Guidelines for community noise. Geneva: World Health Organization; 1999. Accessed 1 June 2024. Available at:[Web]
- [9] Pisani MA, Friese RS, Gehlbach BK et al. Sleep in the intensive care unit. *Am J Respir Crit Care Med*, 2015; 191: 731-738.[DOI]
- [10] Darbyshire JL, Young JD. An investigation of sound levels on intensive care units with reference to the WHO guidelines. *Crit Care*, 2013; 17: R187.[DOI]
- [11] Cordova AC, Logishetty K, Fauerbach J et al. Noise levels in a burn intensive care unit. *Burns*, 2013; 39: 44-48.[DOI]
- [12] Elliott RM, McKinley SM, Eager D. A pilot study of sound levels in an Australian adult general intensive care unit. *Noise Health*, 2010; 12: 26-36.[DOI]
- [13] Morales-Cané I, Moral-Arroyo JA, Debbia F et al. Impact of sound levels on physiological and consciousness state of cardiovascular patients. *Nurs Crit Care*, 2022; 27: 240-250.[DOI]
- [14] Pal J, Taywade M, Pal R et al. Noise pollution in intensive care unit: a hidden enemy affecting the physical and mental health of patients and caregivers. *Noise Health*, 2022; 24: 130-136.[DOI]
- [15] Hu RF, Jiang XY, Chen J et al. Non-pharmacological interventions for sleep promotion in the intensive care unit. *Cochrane Database Syst Rev*, 2015; 10: CD008808.[DOI]
- [16] Gallacher S, Enki D, Stevens S et al. An experimental model to measure the ability of headphones with active noise control to reduce patient's exposure to noise in an intensive care unit. *Intensive Care Med Exp*, 2017; 5: 47.[DOI]
- [17] Xyrichis A, Wynne J, Mackrill J et al. Noise pollution in hospitals. *BMJ*, 2018; 363: k4808.[DOI]
- [18] Konkani A, Oakley B. Noise in hospital intensive care units: a critical review of a critical topic. *J Crit Care*, 2012; 27: e1-e9.[DOI]
- [19] Tegnestedt C, Günther A, Reichard A et al. Levels and sources of sound in the intensive care unit - an observational study of three room types. *Acta Anaesthesiol Scand*, 2013; 57: 1041-1050.[DOI]
- [20] Stafford A, Haverland A, Bridges E. Noise in the ICU. Am J Nurs, 2014; 114: 57-63.[DOI]
- [21] Martinez FE, Poulter AL, Seneviratne C et al. ICU patients' perception of sleep and modifiable versus non-modifiable factors that affect it: a prospective observational study. *J Clin Med*, 2022; 11: 3725.[DOI]
- [22] Mori C, Boss K, Indermuhle P et al. Is it noise? Factors linked with sleep interruption in hospitalized patients. *Clin Nurse Spec*, 2021; 35: 199-207.[DOI]
- [23] Buxton OM, Ellenbogen JM, Wang W et al. Sleep



https://doi.org/10.53964/jmnpr.2024015

disruption due to hospital noises: a prospective evaluation. Ann Intern Med, 2012; 157: 170-179.[DOI]

[24] Gabor JY, Cooper AB, Crombach SA et al. Contribution of

the intensive care unit environment to sleep disruption in mechanically ventilated patients and healthy subjects. Am J Respir Crit Care Med, 2003; 167: 708-715.[DOI]

